



## Density of Pennies

Wear Safety Goggles and an Apron.



### Materials

- balance
- 100 mL graduated cylinder
- 40 pennies dated before 1982
- 40 pennies dated after 1982
- water

### Procedure

1. Using the balance, determine the mass of the 40 pennies minted prior to 1982. Repeat this measurement two more times. Average the results of the three trials to determine the average mass of the pennies.
2. Repeat step 1 with the 40 pennies minted after 1982.
3. Pour about 50 mL of water into the 100 mL graduated cylinder. Record the exact volume of the water. Add the 40 pennies minted before 1982. Record the volume of the water and pennies. Repeat this process two more times. Determine the volume of the pennies for each trial. Average the results of those trials to determine the average volume of the pennies.
4. Repeat step 3 with the 40 pennies minted after 1982.
5. Review your data for any large differences between trials that could increase the error of your results. Repeat those measurements.

6. Use the average volume and average mass to calculate the average density for each group of pennies.
7. Compare the average calculated densities with the density of the copper listed in Table 2-4.

### Discussion

1. Why is it best to use the results of three trials rather than a single trial for determining the density?
2. How did the densities of the two groups of pennies compare? How do you account for any difference?
3. Use the results of this investigation to formulate a hypothesis about the composition of the two groups of pennies. How could you test your hypothesis?

### SAMPLE PROBLEM 2-1

A sample of aluminum metal has a mass of 8.4 g. The volume of the sample is 3.1 cm<sup>3</sup>. Calculate the density of aluminum.

### SOLUTION

Given: mass ( $m$ ) = 8.4 g  
volume ( $V$ ) = 3.1 cm<sup>3</sup>

Unknown: density ( $D$ )

$$\text{density} = \frac{\text{mass}}{\text{volume}} = \frac{8.4 \text{ g}}{3.1 \text{ cm}^3} = 2.7 \text{ g/cm}^3$$